

Packaging MDL Specimens for Local MDL Driver Pick-up.

1

Insert labeled vial into mailer.



Insert the mailer into biohazard bag.

2



3

Be sure to include a completed test requisition form for each specimen in the front pocket of the biohazard bag.



Place specimens in your MDL Specimen Lock box for courier pick-up.

4



Contact us at
877 269 0090 · www.mdlab.com

A MEMBER OF GENESIS BIOTECHNOLOGY GROUP

Packaging MDL Specimens for Local MDL Driver Pick-up.

1
Insert labeled vial
into mailer.



3
Be sure to include a completed
test requisition form for each
specimen in the front pocket of
the biohazard bag.



2
Insert the mailer into
biohazard bag.



4
Place specimens in your
MDL Specimen Lock box for
courier pick-up.



Contact us at 877 269 0090 · www.mdlab.com