

# Packaging MDL Specimens for Local Courier Pick-up.

1 Insert labeled vial into mailer.



2 Insert the mailer into biohazard bag.



3 Be sure to include a completed test requisition form for each specimen in the front pocket of the biohazard bag.



A MEMBER OF GENESIS BIOTECHNOLOGY GROUP

Medical Diagnostic Laboratories  
www.mdlab.com • 877.269.0090



A DIVISION OF

CLINICAL DIAGNOSTICS

MEDICAL DIAGNOSTIC  
LABORATORIES

TOLL FREE 1-877-269-0090

www.mdlab.com

IH0086 Upd: 9\_2023

# Packaging MDL Specimens for Local Courier Pick-up.

Insert labeled vial  
into mailer.



Insert the mailer into  
biohazard bag.



Be sure to include a completed test  
requisition form for each specimen  
in the front pocket of the biohazard  
bag.



Place specimens in your MDL  
Specimen Lock box for courier  
pick-up.



# 4